





SMARTCLIP™
SELF-LIGATING APPLIANCE SYSTEM

HYBRID RECTANGULAR ARCHWIRE



Wire sizes and types are easily identified with hygienic single-use packaging.
Available in Unit of Use packaging only. 10 units per package.

| Dimension | Alloy Type* | U/L | OrthoForm™ I Tapered | OrthoForm™ II Square | OrthoForm™ III Ovoid |
|---------------------------------|-------------|-----|----------------------|----------------------|----------------------|
| .017 x .025 ■ Hybrid | SS | U | 9293-889 | 9293-891 | 9293-893 |
| | SS | L | 9293-890 | 9293-892 | 9293-894 |
| .017 x .025 ■ Hybrid Dimpled | SE | U | 9293-877 | 9293-879 | 9293-881 |
| | SE | L | 9293-878 | 9293-880 | 9293-882 |
| | NC | U | 9293-883 | 9293-885 | 9293-887 |
| | NC | L | 9293-884 | 9293-886 | 9293-888 |
| .018 x .025 ■ Hybrid | SE | U | 9293-903 | 9293-905 | 9293-907 |
| | SE | L | 9293-904 | 9293-906 | 9293-908 |
| | NC | U | 9293-909 | 9293-911 | 9293-913 |
| | NC | L | 9293-910 | 9293-912 | 9293-914 |
| | SS | U | 9293-915 | 9293-917 | 9293-919 |
| | SS | L | 9293-916 | 9293-918 | 9293-920 |
| .019 x .025 ■ Hybrid | NC | U | 9293-951 | 9293-953 | 9293-955 |
| | NC | L | 9293-952 | 9293-954 | 9293-956 |
| | SS | U | 9293-957 | 9293-959 | 9293-961 |
| | SS | L | 9293-958 | 9293-960 | 9293-962 |
| .021 x .025 ■ Hybrid | SE | U | 9293-931 | 9293-933 | 9293-935 |
| | SE | L | 9293-932 | 9293-934 | 9293-936 |
| | SS | U | 9293-943 | 9293-945 | 9293-947 |
| | SS | L | 9293-944 | 9293-946 | 9293-948 |

*NC (Nitinol Classic). SE (Nitinol Super Elastic). SS (Stainless Steel)


3M Science.
Applied to Life.™



SMARTCLIP™
SELF-LIGATING APPLIANCE SYSTEM

HYBRID RECTANGULAR ARCHWIRES

ROUNDED CONTOUR
WIRES FOR ENHANCED
PERFORMANCE



- Use in all instances where standard rectangular wires are used
- Contoured surface area improves sliding mechanics
- Hybrid wire shape facilitates insertion into appliances



3M Oral Care
2510 Conway Avenue
St. Paul, MN 55144-1000 USA
Phone 1-800-423-4588
Web 3M.com/ortho

3M Canada
Health Care Division
300 Tartan Dr.
London, ON N5V 4M9 Canada
Phone 1-800-443-1661

3M, SmartClip, and OrthoForm
are trademarks of 3M.
Used under license in Canada.
© 3M 2016. All rights reserved.
70-2021-3926-8 1601

THE BEST OF BOTH WORLDS IN A HYBRID WIRE


A BETTER COMBINATION

Archwire working force, patient comfort and ease of engagement/disengagement have come together in 3M SmartClip™ Hybrid Rectangular Archwires. By precisely rounding the corners of standard rectangular archwires, SmartClip Hybrid Rectangular Archwires blend the power of a rectangular wire with the reduced binding of a round wire. The hybrid archwire design enhances archwire sliding mechanics and results in improved engagement/disengagement.


FACILITATES WIRE ENGAGEMENT

Unlike the standard rectangular archwires, SmartClip Hybrid Rectangular Archwires present a more radiused corner that is less likely to result in a positive stop when engaging an archwire into the slot. Importantly, the hybrid archwires have a greater tolerance for misalignment of the wire to the bracket slot, and is less likely to encounter contact point interference on engagement.


Wire Insertion into Upper Lateral Brackets



Standard Rectangular Wire




Hybrid Wire



Hybrid Wire


Compared to standard rectangular wire, the shape of SmartClip Hybrid Rectangular Archwire provides more tolerance in the wire alignment to the bracket slot without encountering contact point interference (arrows).



Contact angles (red arrows) create binding

LESS BINDING FOR BETTER SLIDING PERFORMANCE

The rounded edges of the SmartClip Hybrid Rectangular Archwire reduce binding and notching and improve sliding performance without compromising archwire integrity.



Source: 3M lab

With a reduction of the coefficient of binding inherent in SmartClip Hybrid Archwires, retraction mechanics are enhanced and are more effective.


NO NEED TO SACRIFICE TORQUE

SmartClip Hybrid Archwires are offered in a selection of sizes to match your torque requirements, based on optimal slot and archwire interaction. As a result you can achieve similar torque expression with SmartClip Hybrid Rectangular Archwires that you are now achieving with standard rectangular archwires.

A WIRE TO MEET YOUR NEEDS

Truly an archwire to meet all of your needs, the SmartClip Hybrid Rectangular Archwire can be used at any point in the working and finishing stages of treatment. Incorporating the SmartClip Hybrid Rectangular Archwire into your treatment progression can enhance ease of use, reduce binding and can improve patient comfort – a winning combination for your practice and your patients.

WORKS GREAT WITH SMARTCLIP™ SELF-LIGATING APPLIANCES



The Hybrid Archwire offers lower disengagement forces when using SmartClip Self-Ligating Appliances.