

## Transbond<sup>™</sup> Plus Light Cure Band Adhesive Unitek<sup>™</sup> Multi-Cure Glass Ionomer Band Cement

| Characteristic Comparison |  |   |
|---------------------------|--|---|
|                           | Transbond <sup>™</sup> Plus<br>Light Cure Adhesive   | Unitek <sup>™</sup> Multi-Cure Glass Ionomer<br>Band Cement |
|                           | Set Data Set Data   Parale Set Data <t< th=""><th></th></t<> |   |
| Banding Field             |  |   |
| Dry                       | Yes  | Yes   |
| Moist-Water               | No   | Yes   |
| Moist-Saliva              | No   | Yes   |
| Curing                    |  |   |
| Light Cure                | 12 Seconds*<br>(3 seconds per cusp)  | 12 Seconds*<br>(3 seconds per cusp)                         |
| Chemical Cure             | No   | 5 Minutes   |
| Adhesion                  |  |   |
| Mechanical<br>Adhesion    | Yes  | Yes   |
| Chemical Adhesion         | No   | Yes   |
| Delivery                  |  |   |
|                           | Single Paste (No Mixing)   | Powder/Liquid (Mixing Required)                             |
| Storage                   |  |   |
|                           | Refrigerate  | Room Temperature  |

\*Ortholux<sup>™</sup> Luminous Curing Light

| Adhesive Clean-Up        |   |   |  |
|--------------------------|---|---|--|
|                          | On Teeth  | On Instruments  |  |
| Before<br>Adhesive Cures | Use a scaler to remove excess adhesive                        | Clean immediately. Once cured, adhesive can be difficult to remove. |  |
|                          | Use a damp 2×2 to remove excess adhesive off occlusal surface | Remove adhesive with damp cloth or immerse instrument in water      |  |
| After<br>Adhesive Cures  | Dessicate adhesive with air                                   | Abrade adhesive lightly until removed                               |  |
|                          | Use a scaler or ultrasonic scaler                             | Place in ultrasonic cleaner   |  |

## **Band Adhesives**

## Frequently Asked Questions

- 1. What is the purpose of a band adhesive?
- A. The purpose of a band adhesive is to act as a "grout" to fill the irregularities between the band and the tooth surface.
- 2. Why is Unitek<sup>™</sup> Multi-Cure Glass Ionomer Band Cement stronger than Transbond<sup>™</sup> Plus Light Cure Band Adhesive?

A. Unitek Multi-Cure band

to both the band and

cement is stronger because

it provides both chemical

and mechanical adhesion



- enamel surface of the tooth.3. Under which conditions would I need the extra strength of Unitek Multi-Cure GI band cement?
- A. In some offices the extra adhesion offered by a glass ionomer is a welcome advantage where bands fit poorly, where there are isolation issues and in situations where forces greater than those normally exerted from headgear and rapid palatal expanders are expected. The trade-off is increased mixing time, increased clean-up time and increased difficulty when removing the cement from the teeth.
- 4. Under which conditions would I want to use Transbond Plus light cure band adhesive?
- A. In offices where band fit, moisture control and loose bands are not an issue. In these offices, Transbond Plus band adhesive is the best choice for ease-of-use and ease of clean-up.
- 5. What is the effect of moisture on Unitek Multi-Cure GI band cement & Transbond Plus light cure band adhesive?
- A. Hybrid glass ionomers, like Unitek Multi-Cure band cement, work best on a moist surface; however, glass ionomers do not tolerate additional water/saliva introduced <u>during the curing process</u>.

Compomers, like Transbond Plus band adhesive, are not as tolerant of moisture as glass ionomers. An effort should be made to keep the tooth dry when using Transbond Plus band adhesive.

- 6. Will decreasing the time I light cure Transbond Plus band adhesive affect band retention?
- A. Yes, Transbond Plus band adhesive must be cured for the full 12 seconds <u>since there is not a chemical</u> <u>cure component to this adhesive</u>. Archwires may be placed immediately after curing. (Reference time using Ortholux<sup>™</sup> Luminous Curing Light).
- 7. Will decreasing the time I light cure Unitek Multi-Cure GI band cement affect band retention?
- A. No, since the Unitek Multi-Cure band cement also cures chemically, band strength will not be affected if archwires are not placed for 5 minutes. Five minutes is the minimum time needed for the adhesive to fully set-up chemically. If the Unitek Multi-Cure band cement is light cured for the full 12 seconds, archwires can be placed immediately.
- 8. Should I consider using both Transbond Plus light cure band adhesive and Unitek Multi-Cure GI band cement in my orthodontic practice?
- A. Each product provides unique advantages to different banding situations. In some circumstances, one may be more suitable than the other, making both valuable tools in the orthodontic practice.





3M Oral Care 2510 Conway Avenue St. Paul, MN 55144-1000 USA Phone 1-800-423-4588 Web 3M.com/ortho 3M Canada Health Care Division 300 Tartan Dr. London, ON N5V 4M9 Canada Phone 1-800-443-1661

3M, Ortholux, Transbond, and Unitek are trademarks of 3M. Used under license in Canada. © 3M 2016. All rights reserved. 70-2021-3913-8 1601