



Catalog N°7

Elastomerics



Experts in Orthodontics
GC Orthodontics Europe GmbH
www.gc.dental/ortho



Elastomerics

Elastomeric Color Chains.....	191
Orchestra Intra and Extraoral Orchestra elastics Plus	192
Orchestra Intra and Extraoral Elastics.....	193
Japanese Elastomeric Chain – GC Ortho Chain.....	194
GC Ties Elastomeric Ligatures	194
Hollow or Solid Elastomeric Threads.....	195

Elastomeric Color Chains

Available in a wide range of colors that remain stable with time, our color chains are sold in 4.5 meter spools in either long or short.

Color		Short	Long
	Silver	88-1001-0000	88-3001-0000
	Pearl	88-1002-0000	—
	Clear	88-1003-0000	88-3003-0000
	Blue	88-1004-0000	88-3004-0000
	Red	88-1005-0000	—
	Green	88-1006-0000	88-3006-0000
	Black	88-1007-0000	88-3007-0000
	White	88-1008-0000	—
	Pink	88-1009-0000	—
	Metallic Purple	88-1010-0000	—
	Metallic Blue	88-1011-0000	—
	Metallic Red	88-1012-0000	—
	Lavender	88-1013-0000	88-3013-0000
	Fuchsia	88-1014-0000	—
	Yellow	88-1015-0000	—
	Orange	88-1016-0000	88-3016-0000
	Navy Blue	88-1017-0000	—
	Light Blue	88-1018-0000	88-3018-0000
	Maroon	88-1019-0000	—
	Grey	88-1020-0000	88-3020-0000
	Glow Blue	88-1021-0000	—
	Metallic Gold	88-1022-0000	88-3022-0000



Orchestra Intra and Extraoral Orchestra elastics Plus

Our range of intra and extraoral elastics are made from the highest quality medical grade latex.

The stated force corresponds to the force three times its original listed diameter. Intraoral elastics are packaged by 5000 pcs. Extraoral elastics are packaged by 2500 pcs.



Orchestra Intra and Extraoral Elastics

The stated force corresponds to the force three times its original listed diameter. Intraoral elastics are packaged by 5000 pcs. Extraoral elastics are packaged by 2500 pcs. Each bag contain 1 insertion plastic instrument. Latex free elastics are listed on the chart below in green.



EXTRA-ORAL				
Force Size	2.5 oz (70,88 g)	4.5 oz (113,40 g)	6.5 oz (170,10 g)	8,5 oz (240, 98 g)
3 mm (1/8")	85-0321-2000 Harp	85-0331-2000 Xylophone	85-0341-2000 Tuba	
4 mm (3/16")	85-0421-2000 Tambourine	85-0431-2000 Piano	85-0441-2000 Trumpet	85-0461-2000 Trombone
6 mm (1/4")	85-0621-2000 Oboe	85-0631-2000 Flute	85-0641-2000 Baritone	85-0661-2000 French Horn
8 mm (5/16")	85-0821-2000 Cello	85-0831-2000 Clarinet	85-0841-2000 Saxophone	85-0861-2000 Guitar
10 mm (3/8")	85-1021-2000 Cymbals	85-1031-2000 Tympani		85-1061-2000 Violin

EXTRA-ORAL				
Force Size	2.5 oz (70,88 g)	4.5 oz (113,40 g)	6.5 oz (170,10 g)	8,5 oz (240, 98 g)
3 mm (1/8")	85-0321-0000 Harp	85-0331-0000 Xylophone	85-0341-0000 Tuba	
4 mm (3/16")	85-0421-0000 Tambourine	85-0431-0000 Piano	85-0441-0000 Trumpet	85-0461-0000 Trombone
4 mm (3/16") (sans latex)		85-0430-0000 Jazz (latex free)	85-0440-0000 Rock-n-Roll (latex free)	
6 mm (1/4")	85-0621-0000 Oboe	85-0631-0000 Flute	85-0641-0000 Baritone	85-0661-0000 French Horn
6 mm (1/4") (sans latex)		85-0630-0000 Hip-Hop (latex free)	85-0640-0000 Classical (latex free)	
8 mm (5/16")	85-0821-0000 Cello	85-0831-0000 Clarinet	85-0841-0000 Saxophone	85-0861-0000 Guitar
10 mm (3/8")	85-1021-0000 Cymbals	85-1031-0000 Tympani		85-1061-0000 Violin



Japanese Elastomeric Chain – GC Ortho Chain

Made from high-density and high-purity materials, this chain is precision-cut with identical intervals between links. GC Ortho Chain exerts light and continuous force throughout its action, which is essential for optimal tooth movement. Available in 2 m spools in grey or clear.

Color	Short	Middle	Long
Transparent	86-1010-0000	86-1020-0000	86-1030-0000
Grey	86-2010-0000	86-2020-0000	86-2030-0000



GC Ties Elastomeric Ligatures

Latex-free, GC Ties ligatures allow for an upper and lower archwire change with a single stick. The high quality of the polyurethane enables these ligatures to recover 90 % of their initial shape after being stretched.

Available in bags of 50 sticks, 20 ligatures per stick.

Outer diameter: .120"

■	Silver	87-1201-0000	■	Lavender	87-1213-0000
□	Pearl	87-1202-0000	■	Fuchsia	87-1214-0000
□	Clear	87-1203-0000	■	Yellow	87-1215-0000
■	Blue	87-1204-0000	■	Orange	87-1216-0000
■	Red	87-1205-0000	■	Navy Blue	87-1217-0000
■	Green	87-1206-0000	■	Light Blue	87-1218-0000
■	Black	87-1207-0000	■	Maroon	87-1219-0000
□	White	87-1208-0000	■	Grey	87-1220-0000
■	Pink	87-1209-0000	■	Glow Blue	87-1221-0000
■	Metallic Purple	87-1210-0000	■	Golden	87-1222-0000
■	Metallic Blue	87-1211-0000	■	Glow Clear	87-1223-0000
■	Metallic Red	87-1212-0000	■	Tooth	87-1224-0000
			■	Pearl Green	87-1225-0000

Separators



Very strong, latex-free separators with memory effect, these separators are ideal for creating spaces between molars.

Description	P/N
Dark Blue non radio opaque separator, 960 pcs	89-1701-0000
Radio opaque separator, 1000 pcs	89-1502-0000
Radio opaque separator small	89-1503-0000



Hollow or Solid Elastomeric Threads

Our hollow or solid elastomeric threads are available in diameter sizes of .025" and .030" in 15 m spools.

Hollow Threads	P/N
.025" Grey	89-4325-0000
.025" Clear	89-3325-0000

Solid Threads	P/N
.030" Grey	89-4135-0000
.030" Clear	89-3135-0000

