

**3M** Science.  
Applied to Life.™

## Transbond™ Plus Light Cure Band Adhesive



No Mixing.  
No Contamination.  
No Waste.

### **A convenient, single paste, no mix band adhesive from 3M.**

It should come as no surprise that the company that revolutionized orthodontic adhesives with APC™ Adhesive Coated Appliances also brings you a single paste band adhesive. Once again, we've simplified the process. Transbond™ Plus Light Cure Band Adhesive can maximize your productivity by helping you minimize chair time, reduce adhesive waste and by providing superior band retention.



### **Contact 3M today and discover more about Transbond™ Plus Light Cure Band Adhesive.**

It's the smartest call you can make to enhance efficiency and patient satisfaction.

# Make the smart move to Transbond™ Plus Light Cure Band Adhesive.



## A faster, more predictable solution.

As a single paste, the major advantage of Transbond™ Plus Light Cure Band Adhesive is the elimination of mixing, resulting in reduced clean-up, more consistent adhesive performance, and increased office efficiency.

The adhesive kit comes with five convenient syringes and 50 disposable angulated tips. We've eliminated the mixing pads and concerns about cross-contamination, because you apply the adhesive directly onto the band.



## Function you can see.

Unlike other designs, angulated Transbond™ Plus Adhesive disposable tips apply adhesive directly to the band, reducing waste and the chance for cross-contamination.

The distinctive blue color of Transbond Plus Adhesive is another smart feature. It allows you to easily see the adhesive during banding, making clean-up of excess adhesive easy. During de-banding, you can see any adhesive left on the teeth, which facilitates adhesive remnant clean-up, and provides an easier process for the patient. Used with a high performance light cure unit, such as the Ortholux™ Luminous Curing Light, cure time is a mere 12 seconds.



**3M**

3M Oral Care  
2510 Conway Avenue  
St. Paul, MN 55144-1000 USA

Phone 1-800-423-4588  
Web 3M.com/ortho

3M Canada  
Health Care Division  
300 Tartan Dr.  
London, ON N5V 4M9 Canada

Phone 1-800-443-1661

3M, APC, Ortholux, and Transbond are trademarks of 3M.  
Used under license in Canada.  
© 3M 2016. All rights reserved.  
70-2021-3910-4 1601