Transbond[™] XT Light Cure Adhesive

Transbond™ XT Light Cure Adhesive bonds metal and ceramic brackets to tooth surfaces. Available in both syringes and capsules, this adhesive uses light cure adhesive technology to provide you with additional working time to ensure accurate bracket placement.

The viscosity of Transbond XT adhesive was designed to prevent adhesive run-on and bracket skating that saves you money and reduces adhesive waste. 3M's adhesive-filled capsules and dispensing gun also provide for easy application, increased control and convenience. A syringe delivery system is available too. Best of all, the Transbond XT adhesive's quick cure provides for immediate archwire tie-in so your patients and staff can experience shorter bonding appointments.

Transbond[™] XT Light Cure Adhesive Kit in Capsules

25 capsules (0.2g each) 6 ml Transbond XT Primer 1 adhesive dispensing gun 1 brush holder 60 brush inserts 712-030

Transbond XT Adhesive in Capsules

25 capsules (0.2g each) 712-031

Adhesive Dispensing Gun 712-032



Transbond XT Light Cure

Adhesive Kit in Syringes 2 syringes (4g each) 6 ml Transbond XT Primer 1 brush holder 60 brush inserts 712-035

Transbond XT Adhesive in Syringes

4 syringes (4g each) 712-036

Transbond XT Primer Bottle

6 ml 712-034

Transbond XT Etching Gel System (Blue Tint)

(35% Phosphoric Acid) 2 gel etch syringes (3.5g each) 25 disposable dispensing tips 712-039

Gel Etch Syringe Tips 25 disposable tips

712-037 Brush Holder

4/pk 704-046

Brush Tips 60/pk x2 704-047

Command Cure

 Quick ceramic/metal bracket cure

Extended Working Time

 Allows precise bracket placement

Immediate Bond Strength

 Allows archwire to be tied in immediately following cure

Versatile

- Saves time for rebonds
- Efficient bonding of ceramic and metal brackets


Pre-loaded Syringes or Capsules

No waste of materials

Excellent Handling Properties

- No bracket drift
- Easy flash clean-up

Typical Cure Time Necessary Before Archwire Tie-In





Transbond[™] LR Light Cure Adhesive

Our Transbond™ LR Adhesive has been developed specifically for bonding lingual retainers. Easy to apply, it has optimal handling characteristics that allow you to shape and finish the adhesive around lingual retainer wire for maximum patient comfort. Like the Transbond™ XT Adhesive, working time is extended so accurate placement is enhanced and product waste is minimized.

Transbond™ LR Capsule Kit

25 capsules (0.2g each) 6 ml Transbond XT Primer 1 adhesive dispensing gun 1 brush holder 60 brush inserts 712-033

Transbond LR Capsules

25 capsules (0.2g each) 712-038

Adhesive Dispensing Gun

712-032


EZ3 AlastiK™ Holder

Holds lingual retainer wire securely in place during bonding 10/pk 406-044

Transbond™ XT Primer Bottle

6 ml 712-034

Flowable Adhesive, Direct Bond Strength



Source: 3M lab 2005.



Transbond XT Etching Gel System (Blue Tint)

(35% Phosphoric Acid) 2 gel etch syringes (3.5g each) 25 disposable dispensing tips 712-039

Gel Etch Syringe Tips

25 disposable tips 712-037

Brush Holder

4/pk 704-046

Brush Tips

60/pk x2 704-047